


Product Name: Amb a 6

Catalog # AA01-A32U

Storage and Stability: Store at -20°C . To avoid repeated handling and multiple freeze/thaw cycles aliquot product into smaller quantities.

Allergen Name	: Amb a 6
Alternate Names	: LTP, Non-specific lipid-transfer protein, Allergen Amb a VI, Allergen Ra6, Pollen allergen Amb a 6
Source	: <i>Ambrosia artemisiifolia</i>
Accession Numbers	: GenBank: U89793 UniProt: 004004
Host	: <i>E. coli</i>
Formulation	: Protein stored in 1x PBS, 5% trehalose solution
Molecular weight	: Calculated: 10.1 kDa Observed: ~13 kDa
Purity	: The purity was determined to be > 95% by densitometry
Gel Image	180



1. Hiller KM, et al. Cloning and expression of ragweed allergen Amb a 6. *Scand J Immunol.* Jul;48(1):26-36, 1998.
2. Roebber M, et al. Isolation and properties of a new short ragweed pollen allergen, Ra6. *J Immunol.* Aug;131(2):706-11, 1983.